



REF CG01E  
UDI 4262377860014

## Emergency Hemostatic Dressing Instructions for use

### Description

SpeedM® Emergency Hemostatic Dressing is a sterile, 7.5 x 370 cm, radiopaque (detectable with X-rays), Z-fold hemostatic dressing and is packaged for aseptic removal. This product contains Halloysite (mineral).

### Intended Purpose

SpeedM® is an emergency hemostatic dressing intended for temporary treatment of heavily bleeding traumatic wounds. The product is suitable for any patient to whom the intended purpose can be applied. The product is not suitable for children under 6 years of age.

### Indications

SpeedM® Emergency Hemostatic Dressing is intended for temporary treatment of heavily bleeding wounds and traumatic injuries such as cuts, lacerations and abrasions. For use by trained first responders.

### Instructions for use



1.



2.




**Identify injury**, locate active bleedings. Remove the folded gauze from the packaging. **Care for the heaviest bleeding first.**

**Insert** one end of the **gauze down to the deepest point of the wound.** Then stuff and fill the wound cavity tightly from the bottom until the skin surface is reached. Depending on the size of the wound, several dressings can be inserted into the wound cavity. As the dressing is perforated, it does not matter whether the printed or unprinted side comes into contact with the wound first.




3.



4.

### Training video



**Apply direct, firm pressure** for at least **3 minutes** or until bleeding stops.

**Check that the bleeding has stopped.** Then apply a pressure bandage. Seek medical attention immediately and show the package.

**Removing the product:** Carefully remove the dressing from the wound. Rinse the wound with sterile solution.

### Contraindications

SpeedM® is not indicated for use on the eyeball due to the risk of dehydration or pressure damage. The product is not indicated for use in or in direct contact with the central circulatory system. SpeedM® must not be used in patients with a known allergic reaction to kaolins.

### Warnings

The product is non-absorbable and must be removed from the wound prior to wound closure. SpeedM® Emergency Hemostatic Dressing is not intended as a substitute for careful surgical technique, the use of sterile procedures, and the proper use of ligatures or other conventional hemostatic procedures. SpeedM® should be removed from the wound after a maximum of 24 hours.

### Precautions

For external use only. Sterility not guaranteed if packaging is damaged or opened. Dispose of the product if the packaging is damaged. For single use only. Protect from heat. Protect from moisture. Do not resterilize.

### Storage conditions

Keep dry. Protect from heat, including storage in direct sunlight or in direct contact with heat sources. Store at a temperature of 25°C (77°F) or less.

### Sterility and shelf life

The product is sterilized by gamma radiation and is intended for single use only. Do not resterilize. Do not reuse. Do not use if the sterile bag is damaged or opened. Do not use the product after the expiry date indicated on the packaging.

### Notes

If serious incidents or unexpected events occur in connection with SpeedM®, please inform the manufacturer nonconfidentially at [speedcaremineral.de](mailto:speedcaremineral.de) and/or the responsible authority. Ensure that the product is disposed properly in accordance with national regulations.



CE 0482



Manufacturer  
Manufacturing country  
Usable until  
Batch name  
Catalog number  
Radiation sterilized  
Single sterile barrier  
Do not resterilize  
Do not use if the packaging is damaged  
Do not reuse  
Observe instructions for use  
Medical Device  
Unique device identification  
Upper temperature limit

Manufactured by:  
Speed Care Mineral GmbH  
Genzkower Straße 7  
17034 Neubrandenburg  
Germany  
Tel: +49 395 42392-0  
Fax: +49 395 42392-27  
E-Mail: [info@speedcaremineral.de](mailto:info@speedcaremineral.de)  
[www.speedcaremineral.de](http://www.speedcaremineral.de)  
© 2024 Speed Care Mineral GmbH  
All rights reserved.  
SCM2024012 Rev 00 09/2024  
Made in Germany